

# NEOVELLE™

Extra Pure Hyaluronic Acid

## NEOVELLE™ SOFT

• Superficial dermis

Filling of superficial wrinkles and lip enhancement

## NEOVELLE™ DEFINE

• Middle dermis

Redefinition of facial imperfections and filling of medium-deep furrows.

## NEOVELLE™ VOLUME

• Deep dermis

Restoration and creation of facial volumes.

## NEOVELLE™ HYDRO BOOST

• Superficial and middle - deep dermis

Regenerates, volumizes and moisturizes the skin, improving elasticity.



### REGISTERED OFFICE

Via Torconca, 1 - 47842 S. Giovanni in Marignano (RN)

Partita IVA 03189570405

Tel. +39 0541 410286 - Fax + 39 0541954644

info@marllor.com - www.marllor.com

MADE IN ITALY

PREMIUM QUALITY

# NEOVELLE™

Extra Pure Hyaluronic Acid



Pure  
Safe  
Effective...  
What else?

## TREATMENT AREAS

### NEOVELLE™ SOFT

- 1 Cupid bow definition (philtrum)
- 2 Definition lip contours
- 3 Increase lip volume
- 4 Perioral wrinkles

### NEOVELLE™ DEFINE

- 5 Eyebrow lift
- 6 Non surgical rhinoplasty
- 7 Nasolabial fold
- 8 Zygomatic area
- 9 Periorbital area
- 10 Marrionette line

### NEOVELLE™ VOLUME

- 11 Correction of depressed nasal angle
- 12 Volumetric augmentation of cheekbones
- 13 Plastic correction of the chin
- 14 Rhinoplasty

### NEOVELLE™ HYDRO BOOST


- 15 Cheeks
- 16 Periorbital area
- 17 Mandibular area
- 18 Hands
- 19 Décolleté
- 20 Peri-zygomatic



# NEOVELLE™ SOFT



**Composition:** cross-linked hyaluronic acid  
**Hyaluronic Acid Concentration:** 15mg/ml  
**Filling capacity:** ● ○ ○ ○  
**Packaging:** 2 x 1ml  
**Needles:** 2 x 30G 1/2"  
**Estimated duration of results:** up to 6 months



**Injection depth:** superficial dermis




• BEFORE • AFTER

# NEOVELLE™ DEFINE



**Composition:** cross-linked hyaluronic acid  
**Hyaluronic Acid Concentration:** 20mg/ml  
**Filling capacity:** ● ● ○ ○  
**Packaging:** 2 x 1ml  
**Needles:** 2 x 27G 1/2"  
**Estimated duration of results:** up to 8 months



**Injection depth:** middle dermis




• BEFORE • AFTER


# NEOVELLE™ VOLUME



**Composition:** cross-linked hyaluronic acid  
**Hyaluronic Acid Concentration:** 25mg/ml  
**Filling capacity:** ● ● ● ●  
**Packaging:** 2 x 1ml  
**Needles:** 2 x 27G 1/2"  
**Estimated duration of results:** up to 10 months



**Injection depth:** Deep dermis




• BEFORE • AFTER


# NEOVELLE™ HYDRO BOOST



**Composition:** linear Hyaluronic Acid  
**Hyaluronic Acid Concentration:** 20 mg/mL  
**Packaging:** 1 x 2ml  
**Needles:** 2 x 30G 1/2"  
**Estimated duration of results:** more than 2 months



**Injection depth:** superficial and medium-deep dermis



• BEFORE • AFTER