

CLASSIC

Before

After



Nasolabial fold



Non surgical rhinoplasty

By kind permission of Prof. P. Motolese

PLUS

Before

After



Nasolabial furrow

By kind permission of Prof. P. Motolese

LIPS

Before

After



Lip reshaping and volume enhancement

By kind permission of Prof. P. Motolese



Lower third of the face rejuvenation

By kind permission of Prof. P. Motolese

NEOVELLE
Extra Pure Hyaluronic Acid

©PAUCCIMARKETING

Pure
Safe
Effective...
What else?

Regenyal Laboratories Srl - Via Valle Piana, snc
63074 San Benedetto del Tronto (AP) ITALY

 **marllor**
BIOMEDICAL

REGISTERED HEADQUARTERS
Via Don Minzoni, 12 - 47842 S. Giovanni in Marignano (RN) - Partita IVA 03189570405
OPERATIONAL HEADQUARTERS
Via Torconca, 1 - 47842 S. Giovanni in Marignano (RN)
Tel. +39 0541 410286 - Fax +39 0541 954644
info@marllor.com - www.marllor.com


NEOVELLE LINE

Neovelle™ line is composed of high quality and highly purified cross-linked hyaluronic acid.

The residue of BDDE (butanediol diglycidyl ether) is less than 0,1 ppm (0.08). hyaluronic acid Neovelle™ the purest in the world.

The characteristics of the filler Neovelle™ are:

- Not swelling
- High corrective power
- Homogeneous tissue distribution
- Soft and full results



Neovelle™ should be injected in the medium and medium superficial dermis layer without over-correcting (Neovelle Classic)

Ipoderma or over periosteum (Neovelle Plus)

“ Everything you can imagine, nature has already created. ”



NEOVELLE CLASSIC

High purity, manageability and safety. Injected into the dermis medium to medium surface provides the highest power correction and thanks to its peculiar viscosity in tissues also self-limiting minor overcorrection.



NEOVELLE PLUS

Crosslinked and intercalated filler, the residue of BDDE <0.1 PPM. Volumetric and long-lasting.



NEOVELLE LIPS

Soft and adequate viscosity for a homogeneous and uniform distribution in tissues. It gives a balanced and harmonic correction respecting the morphology of the lips.

CLASSIC

- A eyebrow lift
- B non surgical rhinoplasty
- C nasolabial fold
- D zygomatic area
- E periorbital area
- F marrionette line

PLUS

- G correction of depressed nasal angle
- H volumetric augmentation of cheekbones
- I plastic correction of the chin
- L rhinoplasty

LIPS

- M cupid bow definition (philtrum)
- N definition lip contours
- O increase lip volume
- P perioral wrinkles

